



**EN9000** SERIES  
FLOOR LEVEL BED  
BY ALRICK



# EN9000 SERIES FLOOR LEVEL BED



## A bed with endless possibilities

The Alrick EN9000 Extreme Care Series Bed is one of the lowest beds available on the market and is designed for the most vulnerable patients and residents, lowering to floor level to prevent injuries and accidents caused by falls.

All mechanisms and actuators are enclosed beneath the mattress platform, meaning that this is the *first* true level profiling bed, which has *no end pedestals*.

- SWL of 200KG or 300KG
- No End Columns: No Restrictions
- Clean smooth design

### Specification / Features:

- Electric hi-lo, 100mm – 725mm
- Electric 4-section profiling
- Auto-regression backrest for minimising shear forces
- Electric forward and reverse Trendelenburg
- Chair positioning
- Removable head and footboards
- Optional side rails
- Optional central locking
- Overall length 2100mm
- Overall width 900 mm
- Mattress platform dimensions: 1990mm x 900mm
- Mobilisation features
- Unique Dual SWL (200kg / 300kg) feature\*
- A mattress platform free of restrictions

**caremed**  
HEALTHCARE EQUIPMENT

Contact Care-Med to request a demo, trial or more information:

Bericote Park, Bericote Road, Blackdown Leamington Spa, CV32 6QP

Freephone: 08000 191 246 | Email: [info@care-med.co.uk](mailto:info@care-med.co.uk)